

ΩMEGA

SURGICAL INSTRUMENTS



SURGEON PREFERRED TWIST DRILLS

Expect more from your cutting instruments

EXPERIENCE CUTTING EXCELLENCE


OPTIMIZED PERFORMANCE

Omega's precision twist drills are the culmination of decades worth of cutting accessory manufacturing. Our grinding and finishing processes combined with special selections of premium materials allow our drill bits to cut smoothly and efficiently while maintaining a sharp cutting edge.

- ▶ Premium Quality
- ▶ Pinpoint Accuracy
- ▶ Exceptional Cutting Edge
- ▶ Single Pack Sterilized for Immediate Use

TITANIUM NITRIDE

A Titanium Nitride coating is applied to our 2.5mm bits to provide exceptional wear resistance and reduced friction.



WHY
ΩMEGA

FOCUSED SOLELY ON PREMIUM CUTTING ACCESSORIES

It is our belief that every patient deserves the best possible treatment. That is why Omega focuses solely on cutting accessories. We strive to produce the highest quality cutting products available by continually perfecting our line

RIGOROUS QUALITY STANDARDS

Only the highest grades of stainless steel and carbide are selected for Omega's products. Our manufacturing processes, developed over a period of several decades, blend high tech manufacturing with individual craftsmanship. Only instruments passing a rigorous series of inspections receive the Omega symbol identifying them as the industry standard of quality.

DECADES OF MANUFACTURING EXCELLENCE

Established in 1993, Omega Surgical Instruments specializes as a blade and bur cutting accessory company.

- CE Mark Approved
- ISO 13485:2016 Certified
- Manufactured in the USA

REQUEST A QUOTE

To request a quote call **810-695-9800** or email sales@omegasurgical.com

ΩMEGA

SURGICAL INSTRUMENTS

Omega Surgical Instruments
1072 S. Elms Rd., Suite E
Flint, Michigan 48532

Phone

810.695.9800

Email Orders To

omegasurgicalPO@gmail.com

Quote Request

sales@omegasurgical.com



Proudly manufactured
in the United States

Sterile packaged for
immediate use