

Omega Cast Removal Accessories

OMEGA
SURGICAL INSTRUMENTS

Stryker® 986 Cast Saw Filters

Direct Replacement Filters for Stryker® CastVac™


The Omega CC986-500 vacuum filter is a brand-new direct replacement for the Stryker® 986-500 Cast Saw Vacuum filter. It can be used in both the 986 Mobile and 987 Wall Mount vacuums. Replacing the filter every 6 months is recommended based on usage.



Stryker® 940 Cast Saw Blades

Built to last. Omega cast cutter blades are manufactured with special select alloys to ensure maximum durability while remaining cool.

Stryker 840 Cast Saw Blades



The CC1-001 is designed for Plaster Cast removal with the Stryker 840 Cast Cutter



The CC940-25 is designed for plaster cast removal.



The CC2-5X is designed for synthetic Cast removal.



Our Titanium Nitride Coated CC940-35 is designed for harder cast materials



Manufactured in the USA

Contact us



The CC940-23 features a hardening process making it our maximum durability blade.

Omega Surgical Instruments

1072 South Elms Rd. Suite E
Flint, MI 48532

810.695.9800
www.omegasurgical.com

Email Orders to:
po@omegasurgical.com

Bundle & Save


**Large Quantity Discounts
& Recurring Pricing Available.**

Choose Omega today and experience cutting excellence.

- Usually ships in 24 - 48 hours
- Wide selection of products
- Superior Quality

DeSoutter Cast Saw Blades

Our DCC-635-25 cast cutter blades feature a proprietary design to increase airflow and decrease heat generation.



Why Omega

Decades of Manufacturing Excellence

Established in 1993, Omega Surgical Instruments specializes in surgical blade, bur, and twist drill manufacturing.

Focused On Premium Cutting Accessories

It is our belief that every patient deserves the best possible treatment. That is why Omega focuses solely on cutting accessories. We strive to produce the highest quality cutting products available by continually perfecting our line.

Global Leader

Our specialization has allowed us to supply cutting accessories to more than 30 countries worldwide.

Rigorous Quality Standards

Only the highest grades of stainless steel and carbide are selected for Omega's products. Our manufacturing processes, developed over a period of several decades, blend high tech manufacturing with individual craftsmanship. Only instruments passing a rigorous series of inspections receive the Omega symbol identifying them as the industry standard of quality.

Download our Catalog
www.omegasurgical.com/catalog

Omega Surgical Instruments

1072 South Elms Rd. Suite E
Flint, MI 48532

810.695.9800
www.omegasurgical.com

Email Orders to:
po@omegasurgical.com